

## Membrane Screws & Tacks

Membrane Screws and Membrane Tacks are used for the easy and secure attachment of regenerative membranes to the bone during guided bone regeneration procedures. Whether the preference is for fastening membranes by screwing or tapping, Neoss provides a solution for both options.

### Features & benefits

**Flexible** – A solution for soft to hard bone.

**Easy handling** – Fast and easy membrane fastening and removal. For Screws utilise the Neoss Implant Inserter. For Tacks use the Tack positioning tool from Neoss.

**Screw with self-tapping properties** – Minimizes the need for pre-drilling.

**Easy storage** – One Cassette for both Screws and Tacks.




For information on NeoGen<sup>®</sup>, Neoss non-resorbable PTFE membranes, please visit [www.neoss.com](http://www.neoss.com).

### Facts & figures

**Packaging** – 10-pack in a glass vial in a blister. Non-sterile, delivered clean and ready for sterilization. Single use.

**Material** – Titanium alloy grade 5.

**Insertion torque** – 10 Ncm max.



### Assortment

| Art. No. | Description                  |
|----------|------------------------------|
| 64030    | Membrane Tack 3 mm – 10 pcs  |
| 64031    | Tack Positioning Instrument  |
| 64032    | Tack Mallet                  |
| 64034    | Tack and Screw Cassette      |
| 64035    | Membrane Screw 3 mm – 10 pcs |

*All products might not be available in all markets.*

For ordering details on NeoGen<sup>®</sup> membranes and Neoss Implant Inserters please see Neoss Product Catalogue.

See instructions for use for all procedure details.

Neoss Ltd · Windsor House · Cornwall Road · Harrogate · HG1 2PW · United Kingdom · T +44 1423 817733 · F +44 1423 817744 · E [info@neoss.com](mailto:info@neoss.com)  
 Neoss Australia Pty. Ltd · PO Box 404 · New Farm · QLD 4005 · Australia · T +61 7 3216 0165 · F +61 7 3216 0135 · E [info.au@neoss.com](mailto:info.au@neoss.com)