

OVERSIGT

IPD screwdrivers 100% kompatible med +30 implantatsystemer. Se alle online.



Scan QR link til oversigt




Screwdrivers per implant system

Please note that all IPD scan bodies use the same 1.2 Hex and all IPD angulated screws use the TPA screw driver




IPD Scan Abutments (All implant systems)

	18 mm	25 mm	32 mm
1.2 Hex	 PH-20-18	PH-20-25	PH-20-32





IPD Angulated Screw Channel (All implant systems use TPA screwdriver)

	18 mm	25 mm	32 mm
TPA	 KA-CT-18	KA-CT-25	KA-CT-32





Megagen AnyRidge® [WB] (IPD/WB) Ø 4.0

	18 mm	25 mm	32 mm
1.2 Hex	 PH-20-18	PH-20-25	PH-20-32




Camlog® Camlog® [JA] (IPD/JA) Ø 3.3 | Ø 3.8 | Ø 4.3

	18 mm	25 mm	32 mm
1.25 Hex	 PH-25-18	PH-25-25	PH-25-32
TPA	 KA-CT-18	KA-CT-25	KA-CT-32





ASTRA Evolution® [EB] (IPD/EB) Ø 3.6 | Ø 4.2

	18 mm	25 mm	32 mm
1.27 Hex	 PH-27-18	PH-27-25	PH-27-32
TPA	 KA-CT-18	KA-CT-25	KA-CT-32





DENSPLY Xive® Friadent [IA] (IPD/IA) Ø 3.4 | Ø 3.8 | Ø 4.5

	18 mm	25 mm	32 mm
1.2 Hex	 PH-20-18	PH-20-25	PH-20-32
TPA	 KA-CT-18	KA-CT-25	KA-CT-32





ASTRA Osseospeed® [EA] (IPD/EA) Ø 3.0 | Ø 3.5/4.0 | Ø 4.5/5.0

	18 mm	25 mm	32 mm
1.27 Hex	 PH-27-18	PH-27-25	PH-27-32
TPA	 KA-CT-18	KA-CT-25	KA-CT-32





DENSPLY Ankylos® [IB] (IPD/IB) C/X

	18 mm	25 mm	32 mm
1.0 Hex	 PH-00-18	PH-00-25	PH-00-32
TPA	 KA-CT-18	KA-CT-25	KA-CT-32



DIO Dio® [UB] (IPD/UB) NP/RP

	18 mm	25 mm	32 mm
1.2 Hex	 PH-20-18	PH-20-25	PH-20-32
TPA	 KA-CT-18	KA-CT-25	KA-CT-32



Bestilles på
ipd.puredent.dk



OVERSIGT

IPD screwdrivers 100% kompatible med +30 implantatsystemer. Se alle online.



Scan QR
link til oversigt



Nobel Brånemark system® [AA] (IPD/AA) Ø 3.5 | Ø 4.1 | Ø 5.1

18 mm 25 mm 32 mm

1.2 Hex		PH-20-18	PH-20-25	PH-20-32
UNIGRIP		PU-00-18	PU-00-25	PU-00-32
TPA		KA-CT-18	KA-CT-25	KA-CT-32



Straumann Tissue Level® [DA] (IPD/DA) Ø 4.8 RN | Ø 6.5 WN

18 mm 25 mm 32 mm

TORX		PT-00-18	PT-00-25	PT-00-32
TPA		KA-CT-18	KA-CT-25	KA-CT-32



Nobel Active® [AD] (IPD/AD) Ø 3.0 | Ø 3.5 | Ø 4.3

18 mm 25 mm 32 mm

UNIGRIP		PU-00-18	PU-00-25	PU-00-32
TPA		KA-CT-18	KA-CT-25	KA-CT-32



Straumann Bone Level® + BLT® [DB] (IPD/DB) Ø 3.3 NC | Ø 4.1 RC

18 mm 25 mm 32 mm

TORX		PT-00-18	PT-00-25	PT-00-32
TPA		KA-CT-18	KA-CT-25	KA-CT-32



Nobel replace® select [AC] (IPD/AC) Ø 3.5 | Ø 4.3 | Ø 5.0

18 mm 25 mm 32 mm

UNIGRIP		PU-00-18	PU-00-25	PU-00-32
TPA		KA-CT-18	KA-CT-25	KA-CT-32



Straumann BLX® [DC] (IPD/DC) RB/WB

18 mm 25 mm 32 mm

TORX		PT-00-18	PT-00-25	PT-00-32
TPA		KA-CT-18	KA-CT-25	KA-CT-32

Other screwdrivers

20 mm 25 mm 32 mm

Multi-Unit®		LL-AB-20
PSD (Locator®)		LL-PS-00

Screwdriver kits

Get the full overview at the IPD shop:
www.ipd.puredent.dk/screwdrivers/



Bestilles på
ipd.puredent.dk

